

K0284 rabbit polyclonal antibody

Catalog No: #29588

Package Size: #29588-1 50ul #29588-2 100ul


Description

Product Name	K0284 rabbit polyclonal antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB
Species Reactivity	Hu Ms
Specificity	This antibody detects endogenous levels of K0284 at Human/Mouse
Immunogen Description	Synthesized peptide derived from human K0284
Conjugates	Unconjugated
Accession No.	Swiss-Prot#: Q9Y4F5
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20oC/50%1 year

Application Details

WB 1oC/500-2000

Images



Western blot analysis of lysates from MCF-7 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only and is not intended for use in humans or animals.