

OCSTP polyclonal antibody

Catalog No: #29589

Package Size: #29589-1 50ul #29589-2 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	OCSTP polyclonal antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Synthesized peptide derived from human OCSTP
Conjugates	Unconjugated
Concentration	1mg/ml
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000

Images



Western blot analysis of lysates from AD293 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only and is not intended for use in humans or animals.