

CH25H rabbit polyclonal antibody

Catalog No: #29597

Package Size: #29597-1 50ul #29597-2 100ul



Description

Product Name	CH25H rabbit polyclonal antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Synthesized peptide derived from human CH25H
Conjugates	Unconjugated
Accession No.	Swiss-Prot#: O95992
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20oC/±oC/±1 year

Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.