

VWF Polyclonal Antibody

Catalog No: #29598

Package Size: #29598-1 50ul #29598-2 100ul

Description

Product Name	VWF Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	IHC-p IF
Species Reactivity	Hu Ms Rt
Specificity	VWF Polyclonal Antibody detects endogenous levels of protein.
Immunogen Description	Synthesized peptide derived from part region of human protein
Conjugates	Unconjugated
Accession No.	Swiss-Prot#:P04275
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C ± 5°C / 1 year

Application Details

IHC-p 1:50-300

Background

von Willebrand factor(VWF) Homo sapiens This gene encodes a glycoprotein involved in hemostasis. The encoded preproprotein is proteolytically processed following assembly into large multimeric complexes. These complexes function in the adhesion of platelets to sites of vascular injury and the transport of various proteins in the blood. Mutations in this gene result in von Willebrand disease, an inherited bleeding disorder. An unprocessed pseudogene has been found on chromosome 22.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.