

TIMD4 Rabbit Polyclonal Antibody

Catalog No: #29600

Package Size: #29600-1 50ul #29600-2 100ul

Description

Product Name	TIMD4 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Rat
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human TIMD4 (NP_612388.2).
Conjugates	Unconjugated
Other Names	SMUCKLER;TIM4;TIMD4
Accession No.	Swiss_Prot:Q96H15Gene_ID:91937
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20oC±2oC. Avoid freeze / thaw cycles.

Application Details

WB□1:500 - 1:2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.