

GDF-5 Polyclonal Antibody

Catalog No: #29681

Package Size: #29681-1 50ul #29681-2 100ul



SAB[®]

Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description


Product Name	GDF-5 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB,IHC-p,IF
Species Reactivity	Hu Ms Rt
Specificity	GDF-5 Polyclonal Antibody detects endogenous levels of GDF-5 protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from the Internal region of human GDF5.
Conjugates	Unconjugated
Accession No.	Swiss Prot:P43026
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20oC/±oC/±1 year

Application Details


Western Blot: 1/500 - 1/2000

IHC-p: 1/100-1/300

Images



Western Blot analysis of 293 cells using GDF-5 Polyclonal Antibody. Secondary antibody was diluted at 1:20000



Immunohistochemical analysis of paraffin-embedded rat-brain, antibody was diluted at 1:100

Note: This product is for in vitro research use only and is not intended for use in humans or animals.