

SMARCA5 Antibody

Catalog No: #32545

Package Size: #32545-1 50ul #32545-2 100ul



Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	SMARCA5 Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human;Mouse
Specificity	The antibody detects endogenous level of total SMARCA5 protein.
Immunogen Type	Protein
Immunogen Description	Recombinant fusion protein of human SMARCA5
Conjugates	Unconjugated
Target Name	SMARCA5
Other Names	SMARCA5;ISWI;SNF2H;WCRF135;hISWI;hSNF2H
Accession No.	Uniprot:O60264GeneID:8467
SDS-PAGE MW	122 kDa
Formulation	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at +4°C for short term. Store at -20°C for long term. Avoid freeze/thaw cycle.

Application Details

WB 1:500-1:2000;

Images



Western blot analysis of extracts of various cell lines, using SMARCA5 antibody at 1:1000 dilution.

Background

Helicase that possesses intrinsic ATP-dependent nucleosome-remodeling activity.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.