

## C1QA Polyclonal Antibody

Catalog No: #33201

Package Size: #33201-1 50ul #33201-2 100ul


## Description

|                       |                                                                                                                       |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name          | C1QA Polyclonal Antibody                                                                                              |
| Host Species          | Rabbit                                                                                                                |
| Clonality             | Polyclonal                                                                                                            |
| Isotype               | IgG                                                                                                                   |
| Purification          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Applications          | WB;IHC;IF;ELISA                                                                                                       |
| Species Reactivity    | Human;Rat;Mouse                                                                                                       |
| Specificity           | C1QA Polyclonal Antibody detects endogenous levels of protein.                                                        |
| Immunogen Description | Synthesized peptide derived from part region of human protein                                                         |
| Conjugates            | Unconjugated                                                                                                          |
| Accession No.         | Swiss Prot:P02745GenelD:712                                                                                           |
| Calculated MW         | 26kD                                                                                                                  |
| Concentration         | 1 mg/ml                                                                                                               |
| Formulation           | Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.                                                        |
| Storage               | Store at -20°C/1 year                                                                                                 |


## Application Details

WB 1:500-2000 IHC-p 1:50-200 ELISA 1:5000-20000

## Images



Western blot analysis of lysates from DU145 cells, primary antibody was diluted at 1:1000, 4°over night



Immunohistochemical analysis of paraffin-embedded human tonsil.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.