

CD3 (ABT-CD3) mouse mAb

Catalog No: #33204

Package Size: #33204



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	CD3 (ABT-CD3) mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, Kappa
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Applications	IHC;IF
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of human CD3. Heat-induced epitope retrieval (HIER) Citrate buffer of pH6.0 was highly recommended as antigen repair method in paraffin section
Immunogen Description	Synthesized peptide derived from human CD3
Conjugates	Unconjugated
Target Name	CD3E T3E
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

IHC-p 1:100-500oO WB 1:500-2000. IF 1:50-200

Images



Western blot analysis of CD3Antibody at 1:1000 dilution.



Immunohistochemical analysis of paraffin-embedded
Spleen-high magnification.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.