

EBFP Rabbit Polyclonal Antibody

Catalog No: #38087

Package Size: #38087 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	EBFP Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Applications	WB
Specificity	The EBFP Rabbit Polyclonal Antibody detects EBFP & EBFP Tag protein.
Immunogen Description	Recombinant Protein of BFP
Conjugates	Unconjugated
Target Name	EBFP
Other Names	BFP; EBFP; Blue Fluorescent Protein; Enhanced Blue Fluorescent Protein
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C

Application Details

WB 1:5000;

Images



Western blot analysis of EBFP recombinant protein, using #38087 diluted at 1:5,000.

Background

Enhanced Blue Fluorescent Protein (EBFP) emits a strong blue fluorescence. EBFP can be used as protein “tags” to study the sub-cellular localization of proteins and/or their translocation upon stimulation or as markers for transfection in transient and stable expression systems.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.