

Flotillin-1 Mouse Monoclonal Antibody(6C10)

Catalog No: #38095

Package Size: #38095 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Flotillin-1 Mouse Monoclonal Antibody(6C10)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	6C10
Purification	Affinity purification using immunogen.
Applications	WB
Species Reactivity	Rt Ms
Specificity	Antibody detects endogenous Flotillin-1 protein.
Conjugates	Unconjugated
Target Name	Flotillin-1
Other Names	FLOT 1 antibody, Integral membrane component of caveolae antibody, Reggie 2 antibody
Accession No.	Swiss-Prot#:O75955
SDS-PAGE MW	49kd
Concentration	1.0mg/ml
Formulation	Mouse IgG1 in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000~1:3000

Images



Western blot analysis of 1) Mouse Brain Tissue , 2) Rat Brain tissue, using #38095 diluted at 1:2,000.

Background

Two flotillin members have been described, flotillin-1 and flotillin-2. May act as a scaffolding protein within caveolar membranes, functionally participating in formation of caveolae or caveolae-like vesicles.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.