

AMBRA1 antibody

Catalog No: #38182

Package Size: #38182-1 50ul #38182-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	AMBRA1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total AMBRA1 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human Ambra1 (NP_060219.2).
Conjugates	Unconjugated
Target Name	AMBRA1
Other Names	AMBRA1;DCAF3;WDR94;FLJ20294
Accession No.	Uniprot:Q9C0C7GeneID:55626
SDS-PAGE MW	143kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB□1:500 - 1:2000IHC□1:50 - 1:200


Images



Immunohistochemistry of paraffin-embedded rat lung using AMBRA1 antibody.



Western blot analysis of extracts of various cell lines, using Ambra1 Antibody.



Immunohistochemistry of paraffin-embedded mouse lung using AMBRA1 antibody.

Published Papers

Ye Liuqi;Lin Danlei;Zhang Wen;Chen Shiji;Zhen Yumiao;Akkouche Sara;Liang Xiaoxu;Chong Cheong-Meng;Zhong Hai-Jing; et al., AMBRA1 drives gastric cancer progression through regulation of tumor plasticity, , (2024)

PMID:

Note: This product is for in vitro research use only and is not intended for use in humans or animals.