

ENKD1 antibody

Catalog No: #40407

Package Size: #40407-1 50ul #40407-2 100ul


Description

Product Name	ENKD1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	WB, IHC
Species Reactivity	Hu, Ms
Immunogen Description	Fusion protein of human ENKD1
Conjugates	Unconjugated
Target Name	ENKD1
Other Names	C16orf48; DAKV6410
Accession No.	NCBI Protein#:BC008284
Calculated MW	39 kDa
Concentration	1.4mg/ml
Formulation	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at -20°C/1 year


Application Details

WB dilution:1:500-1:2000IHC dilution:1:50-1:200

Images



Gel: 8%SDS-PAGE, Lysate: 40 µg, Lane 1-3: Human cerebrum tissue, Human cerebella tissue, Mouse brain tissue lysates, Primary antibody:ENKD1 Antibody at dilution 1/1150, Secondary antibody: Goat anti rabbit IgG at 1/5000 dilution, Exposure time: 3 minutes



The image is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using (ENKbody) at dilution 1/80.

Background

The exact function of ENKD1 remains unknown.

Published Papers

Tianyu Zhang;Yixuan Wang;Xiaotong Nie;Xiangrong Chen;Yueyi Jin;Lulu Sun;Ruqian Yang;Jie Wang;Wenqing Xu;Ting Song;Wei Xie;Xiangfeng Chen;Chaojun Li;Jun Zhou;Sijin Wu;Yan Li;Tianliang Li et al., ENKD1 modulates innate immune responses through enhanced geranylgeranyl pyrophosphate synthase activity., , (2025)

[PMID:40048432](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.