

Kif 7 Monoclonal Antibody

Catalog No: #40437

Package Size: #40437-1 50ul #40437-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Kif 7 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	MR4515
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous Kif 7 proteins.
Conjugates	Unconjugated
Target Name	Kif 7
Accession No.	Swiss-Prot: Q2M1P5NCBI Gene ID: 374654
SDS-PAGE MW	360kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Immunohistochemistry: 1:200

Images



IHC staining of Mouse Kidney tissue, diluted at 1:200.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.