

CD85g Polyclonal Antibody

Catalog No: #41665

Package Size: #41665-1 50ul #41665-2 100ul

Description

Product Name	CD85g Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA IHC IF
Species Reactivity	Human Mouse Rat
Specificity	CD85g Polyclonal Antibody detects endogenous levels of CD85g protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human CD85g.
Conjugates	Unconjugated
Target Name	CD85g
Other Names	LILRA4; ILT7; Leukocyte immunoglobulin-like receptor subfamily A member 4; CD85 antigen-like family member G; Immunoglobulin-like transcript 7; ILT-7; CD85g
Accession No.	Swiss-Prot: P59901NCBI Gene ID: 23547
SDS-PAGE MW	55kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

IHC-p;IF(paraffin section) IHC-p: 1:100-300


ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of MCF7 cells using CD85g Polyclonal Antibody



Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:100

Published Papers

Alexander S. P. H.;Wen Jianpei;Li Huihui;Zhou Yufei;Du Hengzhi;Hu Guo;Wen Zheng;Tang Du;Wang Yanwen;Cui Xinwu;Zhou Zhou;Wang Dao Wen;Chen Chen et al., Immunoglobulin attenuates fulminant myocarditis by inhibiting overactivated innate immune response, , (2024)

PMID:

Note: This product is for in vitro research use only and is not intended for use in humans or animals.