

CDO1 Polyclonal Antibody

Catalog No: #44010

Package Size: #44010-1 50ul #44010-2 100ul



SAB[®]

Signalway Antibody

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	CDO1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB
Species Reactivity	Human,Rat,Mouse
Specificity	CDO1 Polyclonal Antibody detects endogenous levels of protein.
Immunogen Description	Synthesized peptide derived from part region of human protein
Conjugates	Unconjugated
Accession No.	SwissProt:Q16878
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20oC/±oC/±1 year

Application Details

WB 1:500-2000

Images



Western blot analysis of lysates from DU145 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only and is not intended for use in humans or animals.