

PSMD12 Polyclonal Antibody

Catalog No: #44017

Package Size: #44017-1 50ul #44017-2 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description


Product Name	PSMD12 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Specificity	PSMD12 Polyclonal Antibody detects endogenous levels of PSMD12 protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human PSMD12.
Conjugates	Unconjugated
Accession No.	Swiss-Prot#: O00232
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000

Immunohistochemistry: 1/100 - 1/300

Images



Western Blot analysis of HepG2 cells using PSMD12 Polyclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.