

GFP-Tag Monoclonal Antibody(2G6)

Catalog No: #44021

Package Size: #44021-1 50ul #44021-2 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	GFP-Tag Monoclonal Antibody(2G6)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	2G6
Isotype	IgG1
Applications	WB
Specificity	The antibody detects GFP, RFP, GFP&RFP fusion proteins.
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Conjugates	Unconjugated
Target Name	GFP-Tag
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:500-1:2000

Images



2ug GFP fusion protein+ Primary antibody dilution at
1) 1:5000
2) 1:10000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.