

Myc-Tag Monoclonal Antibody(5D6)

Catalog No: #44026

Package Size: #44026-1 50ul #44026-2 100ul



SAB[®]

Signalway Antibody

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Myc-Tag Monoclonal Antibody(5D6)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	5D6
Isotype	IgG1
Applications	WB
Specificity	Antibody can detect C-Myc fusion proteins.
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Conjugates	Unconjugated
Target Name	Myc-Tag
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:500-1:2000

Images



Note: This product is for in vitro research use only and is not intended for use in humans or animals.