

Flag-Tag Monoclonal Antibody(2C5)

Catalog No: #44027

Package Size: #44027-1 50ul #44027-2 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Flag-Tag Monoclonal Antibody(2C5)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	2C5
Isotype	IgG1
Applications	WB
Specificity	The antibody detects C-terminal, internal, and N-terminal Flag-tag fusion proteins
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Conjugates	Unconjugated
Target Name	Flag-Tag
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:500-1:2000

Images



Western blot analysis of Recombinant protein, diluted at
1) 1:5000
2) 1:10000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.