

Plant actin Monoclonal Antibody(Q30)

Catalog No: #44034

Package Size: #44034-1 50ul #44034-2 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Plant actin Monoclonal Antibody(Q30)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	Q30
Isotype	IgG1
Applications	WB
Species Reactivity	Plant
Specificity	The antibody detects endogenous plant actin protein.
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Conjugates	Unconjugated
Target Name	Plant actin
Calculated MW	43kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:500-1:2000

Images



Western blot analysis of Arabidopsis, (Q30) diluted at 1:5000

Published Papers

el at., Exogenous Salicylic Acid Regulates Fruiting Body Development, Secondary Metabolite Accumulation, Cell Wall Integrity, and Endogenous Salicylic Acid Content under Heat Stress in Pleurotus ostreatus. In J Agric Food Chem on 2024 Nov 13 by Yanru Hu, Haolan Chen,et al..PMID:39496504, , (2024)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.