

Plant actin Monoclonal Antibody(2E3)

Catalog No: #44045

Package Size: #44045-1 50ul #44045-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Plant actin Monoclonal Antibody(2E3)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	2.00E+03
Isotype	IgG1
Applications	WB
Species Reactivity	Plant
Specificity	The antibody detects endogenous various plants protein.
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Conjugates	Unconjugated
Target Name	Plant actin
Calculated MW	43kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:500-1:2000

Images



Western blot analysis of Arabidopsis, diluted at

1) 1:5000

2) 1:10000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.