

ABCG5 Polyclonal Antibody

Catalog No: #44048

Package Size: #44048-1 50ul #44048-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ABCG5 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB;IF/ICC;ELISA
Species Reactivity	Human,Mouse,Rat
Specificity	ABCG5 Polyclonal Antibody detects endogenous levels of protein.
Immunogen Description	Synthesized peptide derived from human protein .
Conjugates	Unconjugated
Accession No.	Swiss-Prot#: Q9H222
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/5000.

IF 1:100-300

Note: This product is for in vitro research use only and is not intended for use in humans or animals.