

# mStrawberry Monoclonal Antibody(Mix)

Catalog No: #44061

Package Size: #44061-1 50ul #44061-2 100ul


## Description

|                       |                                                                                                            |
|-----------------------|------------------------------------------------------------------------------------------------------------|
| Product Name          | mStrawberry Monoclonal Antibody(Mix)                                                                       |
| Host Species          | Mouse                                                                                                      |
| Clonality             | Monoclonal                                                                                                 |
| Clone No.             | Mix                                                                                                        |
| Isotype               | IgG1                                                                                                       |
| Applications          | WB                                                                                                         |
| Specificity           | The antibody detects mStrawberry & mStrawberry Tag protein. It May had cross reactivity with mRFP protein. |
| Immunogen Type        | Peptide                                                                                                    |
| Immunogen Description | Synthetic Peptide                                                                                          |
| Conjugates            | Unconjugated                                                                                               |
| Target Name           | mStrawberry                                                                                                |
| Concentration         | 1mg/ml                                                                                                     |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                    |
| Storage               | Store at -20°C/1 year                                                                                      |

## Application Details

Western blotting: 1:500-1:2000

## Images



Western blot analysis of mStrawberry recombinant protein, diluted at 1:5000.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.