

β III tubulin Monoclonal Antibody(Mix)

Catalog No: #44062

Package Size: #44062-1 50ul #44062-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	β III tubulin Monoclonal Antibody(Mix)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	Mix
Isotype	IgG1
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous β III tubulin protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Conjugates	Unconjugated
Target Name	β III tubulin
Other Names	TUBB3; TUBB4; Tubulin beta-3 chain; Tubulin beta-4 chain; Tubulin beta-III
Accession No.	Swiss-Prot:Q13509Gene ID:10381
Calculated MW	50kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:500-1:2000

Images



Western blot analysis of 1) Hela, 2) Mouse Brain, 3) Rat Brain, diluted at 1:5000.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.