

## EYFP Monoclonal Antibody(Mix)

Catalog No: #44064

Package Size: #44064-1 50ul #44064-2 100ul


## Description

|                       |                                                                         |
|-----------------------|-------------------------------------------------------------------------|
| Product Name          | EYFP Monoclonal Antibody(Mix)                                           |
| Host Species          | Mouse                                                                   |
| Clonality             | Monoclonal                                                              |
| Clone No.             | Mix                                                                     |
| Isotype               | IgG1                                                                    |
| Applications          | WB                                                                      |
| Specificity           | The antibody detects EYFP and YFP tag fusion proteins.                  |
| Immunogen Type        | Recombinant protein                                                     |
| Immunogen Description | Recombinant Protein                                                     |
| Conjugates            | Unconjugated                                                            |
| Target Name           | EYFP                                                                    |
| Concentration         | 1mg/ml                                                                  |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage               | Store at -20°C/1 year                                                   |

## Application Details

Western blotting: 1:500-1:2000

## Images



Western blot analysis of recombinant EYFP protein, diluted at 1:10000.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.