

PTPRF Polyclonal Antibody

Catalog No: #44068

Package Size: #44068-1 50ul #44068-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PTPRF Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB
Species Reactivity	Human,Mouse,Rat
Specificity	PTPRF Polyclonal Antibody detects endogenous levels of protein.
Immunogen Description	Synthesized peptide derived from part region of human protein
Conjugates	Unconjugated
Accession No.	Swiss-Prot#: P10586 NCBI Gene ID:5792
Calculated MW	209
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000 ELISA 1:5000-20000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.