

HSC70 Monoclonal antibody(6C7)

Catalog No: #44077

Package Size: #44077-1 50ul #44077-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	HSC70 Monoclonal antibody(6C7)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	6C7
Isotype	IgG1
Applications	WB IHC
Species Reactivity	Hu Rt Ms
Specificity	The antibody detects endogenous HSC70 protein
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Conjugates	Unconjugated
Target Name	HSC70
Other Names	HSC70
Accession No.	Swiss-Prot:P11142Gene ID:3312
Calculated MW	70kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western blotting: 1:500-1:2000

Immunohistochemistry: 1:100-1:300

Images



Western blot analysis of 1) HeLa, 2) Mouse Brain Tissue, 3) Rat Brain Tissue with HSC70 Mouse mAb diluted at 1:2,000.



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma using HSC70 Mouse mAb diluted at 1:200.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.