

HSP90 α Monoclonal antibody(7D6)

Catalog No: #44078

Package Size: #44078-1 50ul #44078-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	HSP90 α Monoclonal antibody(7D6)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	7D6
Isotype	IgG1
Applications	WB IHC
Species Reactivity	Hu Rt Ms
Specificity	The antibody detects endogenous HSP90 α protein
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Conjugates	Unconjugated
Target Name	HSP90 α
Other Names	HSP90 α
Accession No.	Swiss-Prot:P07900Gene ID:3320
Calculated MW	95kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western blotting: 1:500-1:2000

Immunohistochemistry: 1:100-1:300

Images



Western blot analysis of 1)Hela, 2)Mouse Brain Tissue, 3) Rat Brain Tissue with HSP90 α Mouse mAb diluted at 1:2,000.



Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma using HSP90 α Mouse mAb diluted at 1:200.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.