

MEK2 Monoclonal Antibody(7D11)

Catalog No: #44083

Package Size: #44083-1 50ul #44083-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MEK2 Monoclonal Antibody(7D11)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	7D11
Isotype	IgG1
Applications	WB
Species Reactivity	Hu Rt Ms
Specificity	MEK2 Mouse Monoclonal Antibody(7D11) detects endogenous levels of MEK2
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Conjugates	Unconjugated
Target Name	MEK2
Accession No.	Swiss-Prot:P36507Gene ID:5605
Calculated MW	45kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:500-1:2000

Images



Western blot analysis of 1) Hela, 2)3T3, 3) Rat Brain Tissue with MEK2 Mouse mAb diluted at 1:2,000.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.