

Rubisco (Large Chain) Monoclonal Antibody(1H1)

Catalog No: #44092

Package Size: #44092-1 50ul #44092-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Rubisco (Large Chain) Monoclonal Antibody(1H1)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	1H1
Isotype	IgG1
Applications	WB
Species Reactivity	Plant
Specificity	Rubisco(Large Chain) Mouse Monoclonal Antibody(1H1) detects endogenous levels of Rubisco(Large Chain)
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Conjugates	Unconjugated
Target Name	Rubisco (Large Chain)
Calculated MW	53kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:2000-1:5000

Images



Western blot analysis of Arabidopsis with Rubisco(Large Chain) Mouse mAb diluted at

1)1:2,000

2)1:5,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.