

Desmin Monoclonal Antibody(1B12)

Catalog No: #44095

Package Size: #44095-1 50ul #44095-2 100ul


Description

Product Name	Desmin Monoclonal Antibody(1B12)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	1B12
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Rt Ms
Specificity	Desmin Mouse Monoclonal Antibody(1B12) detects endogenous levels of Desmin
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Conjugates	Unconjugated
Target Name	Desmin
Accession No.	Swiss-Prot:P17661Gene ID:1674
Calculated MW	53kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Mouse Skeletal Muscle Tissue using DesminBEB™mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Mouse Heart Tissue using DesminMouse mAb diluted at 1:200.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.