

IκB β Monoclonal Antibody(1F3)

Catalog No: #44103

Package Size: #44103-1 50ul #44103-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	IκB β Monoclonal Antibody(1F3)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	1F3
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Rt Ms
Specificity	IκB β Mouse Monoclonal Antibody(1F3) detects endogenous levels of IκB β
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Conjugates	Unconjugated
Target Name	IκB β
Accession No.	Swiss-Prot:Q15653Gene ID:4793
Calculated MW	48kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma using IκB β(mAb diluted at 1:200).



Immunohistochemical analysis of paraffin-embedded Human Kidney Tissue using IκB β Mouse mAb diluted at 1:200.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.