

Transferrin Monoclonal Antibody(5A4)

Catalog No: #44108

Package Size: #44108-1 50ul #44108-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Transferrin Monoclonal Antibody(5A4)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	5A4
Isotype	IgG1
Applications	WB
Species Reactivity	Hu
Specificity	Transferrin Mouse Monoclonal AntibodyB£B`5A4B£B© detects endogenous levels of Transferrin
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Conjugates	Unconjugated
Target Name	Transferrin
Accession No.	Swiss-Prot:P02787Gene ID:7018
Calculated MW	77kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:2000-1:5000

Images



Western blot analysis of Human Serum with Transferrin
Mouse mAb diluted at 1:2,000.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.