

FBN1 Polyclonal Antibody

Catalog No: #44111

Package Size: #44111-1 50ul #44111-2 100ul

Description

Product Name	FBN1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC-p;IF
Species Reactivity	Human;Mouse;Rat
Immunogen Description	The antiserum was produced against synthesized peptide derived from human Fibrillin-1.
Conjugates	Unconjugated
Other Names	FBN1; FBN; Fibrillin-1
Accession No.	Swiss-Prot#: P35555 NCBI Gene ID:2200
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Immunohistochemistry: 1/100 - 1/300

Note: This product is for in vitro research use only and is not intended for use in humans or animals.