

Transferrin Monoclonal Antibody(5F9)

Catalog No: #44117

Package Size: #44117-1 50ul #44117-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Transferrin Monoclonal Antibody(5F9)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	5F9
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu
Specificity	Transferrin Mouse Monoclonal AntibodyB£B"5F9B£B© detects endogenous levels of Transferrin
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Conjugates	Unconjugated
Target Name	Transferrin
Accession No.	Swiss-Prot:P02787Gene ID:7018
Calculated MW	77kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded human-liver using antibody diluted at 1:50.



IHC: 1:100-1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.