

Cleaved Caspase-7 (Asp198) Antibody

Catalog No: #44129

Package Size: #44129-1 50ul #44129-2 100ul

Description


Product Name	Cleaved Caspase-7 (Asp198) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB;IHC-p;IF
Species Reactivity	Human;Mouse
Immunogen Description	The antiserum was produced against synthesized peptide derived from human Caspase 7.
Conjugates	Unconjugated
Accession No.	Swiss Prot:P55210 GeneID:840
Calculated MW	20KD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000

Immunohistochemistry: 1/100 - 1/300

Images



Western Blot analysis of various cells using Cleaved Caspase-7 (Asp198) Antibody



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using Cleaved Caspase-7 (Asp198) Antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.