

P70 S6 Kinase Mouse Monoclonal Antibody(1C7)

Catalog No: #44131

Package Size: #44131-1 50ul #44131-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	P70 S6 Kinase Mouse Monoclonal Antibody(1C7)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	1C7
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Rt Ms
Specificity	P70 S6 Kinase Mouse Monoclonal Antibody(1C7) detects endogenous levels of RPS6KB1
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Conjugates	Unconjugated
Target Name	P70 S6 Kinase
Other Names	RPS6KB1
Accession No.	Swiss-Prot:P23443Gene ID:6198
Calculated MW	70,85kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma Tissue using P70 S6 Kinase Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using P70 S6 Kinase Mouse mAb diluted at 1:200.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.