

STAT1 Mouse Monoclonal Antibody(8H11)

Catalog No: #44134

Package Size: #44134-1 50ul #44134-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	STAT1 Mouse Monoclonal Antibody(8H11)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	8H11
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Rt Ms
Specificity	STAT1 Mouse Monoclonal Antibody(8H11) detects endogenous levels of STAT1
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Conjugates	Unconjugated
Target Name	STAT1
Other Names	STAT1
Accession No.	Swiss-Prot:P42224Gene ID:6772
Calculated MW	84,91kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Colon Tissue using STAT1 Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using STAT1 Mouse mAb diluted at 1:200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.