

JAK1 Mouse Monoclonal Antibody(7G6)

Catalog No: #44135

Package Size: #44135-1 50ul #44135-2 100ul


Description

Product Name	JAK1 Mouse Monoclonal Antibody(7G6)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	7G6
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Rt Ms
Specificity	JAK1 Mouse Monoclonal Antibody(7G6) detects endogenous levels of JAK1
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Conjugates	Unconjugated
Target Name	JAK1
Other Names	JAK1
Accession No.	Swiss-Prot:P23458Gene ID:3716
Calculated MW	130kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma Tissue using JAK1 Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using JAK1 Mouse mAb diluted at 1:200.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.