

α-tubulin (Acetyl Lys40) Monoclonal Antibody

Catalog No: #44147

Package Size: #44147-1 50ul #44147-2 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	α-tubulin (Acetyl Lys40) Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Applications	WB;IF/ICC;IHC-p
Species Reactivity	Human;Rat;Mouse;Drosophila
Immunogen Description	Synthetic Peptide of α-tubulin (Acetyl Lys40)
Conjugates	Unconjugated
Accession No.	Swiss Prot:P68363GeneID:7277
Calculated MW	52KD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

WB 1:1000-2000


IF 1:200

IHC 1:50-100


Images



Western blot analysis of extracts from Hela cells



Immunohistochemical analysis of paraffin-embedded Human-uterus tissue.



Immunofluorescence analysis of Human-liver-cancer tissue.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.