

GSDME n-ternal rabbit pAb

Catalog No: #44148

Package Size: #44148-1 50ul #44148-2 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	GSDME n-ternal rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC-p;IF/ICC
Species Reactivity	Human;Mouse;Rat
Specificity	The antibody detects endogenous full lenth and n-ternal fragment of gsdme protein,
Immunogen Description	Synthesized peptide derived from human GSDME.
Conjugates	Unconjugated
Accession No.	Swiss Prot:O60443
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

WB: 1:5000-20000


IHC: 1:200-300

IF: 1:50-200

Images



Western blot analysis of lysates from 1)mouse-brain, 2)Hela cells treated by UV 15min,3) HepG2 cells



Immunofluorescence analysis of un-treated (left) A549 and UV treated(right) A549 cell.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.