

Cystatin C Mouse Monoclonal Antibody(3B12)

Catalog No: #44158

Package Size: #44158-1 50ul #44158-2 100ul


Description

Product Name	Cystatin C Mouse Monoclonal Antibody(3B12)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	3B12
Isotype	IgG1
Applications	WB
Species Reactivity	Hu
Specificity	Cystatin C Mouse Monoclonal Antibody(3B12) detects endogenous levels of CST3
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Conjugates	Unconjugated
Target Name	Cystatin C
Other Names	CST3
Accession No.	Swiss-Prot:P01034Gene ID:1471
Calculated MW	14kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:1000-1:2000

Images



Western Blot analysis of Cystatin C protein using antibody diluted at 1:1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.