

Gamma Tubulin Mouse Monoclonal Antibody(6C12)

Catalog No: #44180

Package Size: #44180-1 50ul #44180-2 100ul

Description


Product Name	Gamma Tubulin Mouse Monoclonal Antibody(6C12)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	6C12
Isotype	IgG1
Applications	WB IHC
Species Reactivity	Hu Rt Ms
Specificity	Gamma Tubulin Mouse Monoclonal Antibody(6C12) detects endogenous levels of TUBG1
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Conjugates	Unconjugated
Target Name	Gamma Tubulin
Other Names	TUBG1
Accession No.	Swiss-Prot:P23258Gene ID:7283
Calculated MW	50kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western blotting: 1:1000-1:2000


Immunohistochemistry: 1:100-1:200

Images





Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma Tissue using Gamma Tubulin Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using Gamma Tubulin Mouse mAb diluted at 1:200.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.