

Luciferase Mouse Monoclonal Antibody(9D6)

Catalog No: #44185

Package Size: #44185-1 50ul #44185-2 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Luciferase Mouse Monoclonal Antibody(9D6)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	9D6
Isotype	IgG1
Applications	WB
Species Reactivity	Dm
Specificity	Luciferase Mouse Monoclonal Antibody(9D6) detects endogenous levels of Luciferase
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Conjugates	Unconjugated
Target Name	Luciferase
Other Names	Luciferase
Accession No.	Swiss-Prot:P08659Gene ID:
Calculated MW	60kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:1000-1:2000

Images



Western blot analysis of Firefly Luciferase Transfected Hela Cell Lysate using Luciferase Mouse mAb diluted at 1:2000.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.