

PTEN Mouse Monoclonal Antibody(2C10)

Catalog No: #44186

Package Size: #44186-1 50ul #44186-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PTEN Mouse Monoclonal Antibody(2C10)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	2C10
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu
Specificity	PTEN Mouse Monoclonal Antibody(2C10) detects endogenous levels of PTEN
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Conjugates	Unconjugated
Target Name	PTEN
Other Names	PTEN
Accession No.	Swiss-Prot:P60484Gene ID:5728
Calculated MW	54kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma Tissue using PTEN Mouse mAb diluted at 1:200.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.