

Epsilon Tubulin Mouse Monoclonal Antibody(3C1)

Catalog No: #44189

Package Size: #44189-1 50ul #44189-2 100ul


Description

Product Name	Epsilon Tubulin Mouse Monoclonal Antibody(3C1)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	3C1
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Rt Ms
Specificity	Epsilon Tubulin Mouse Monoclonal Antibody(3C1) detects endogenous levels of Epsilon Tubulin
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Conjugates	Unconjugated
Target Name	Epsilon Tubulin
Other Names	Epsilon Tubulin
Accession No.	Swiss-Prot:Q9UJT0Gene ID:51175
Calculated MW	55kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma Tissue using Epsilon Tubulin Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using Epsilon Tubulin Mouse mAb diluted at 1:200.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.