

MMP2 Mouse Monoclonal Antibody(1H1)

Catalog No: #44194

Package Size: #44194-1 50ul #44194-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MMP2 Mouse Monoclonal Antibody(1H1)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	1H1
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu Rt Ms
Specificity	MMP2 Mouse Monoclonal Antibody(1H1) detects endogenous levels of MMP2
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Conjugates	Unconjugated
Target Name	MMP2
Other Names	MMP2
Accession No.	Swiss-Prot:P08253Gene ID:4313
Calculated MW	64B£B~72kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma Tissue using MMP2 Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using MMP2 Mouse mAb diluted at 1:200.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.