

FILA Polyclonal Antibody

Catalog No: #44198

Package Size: #44198-1 50ul #44198-2 100ul


Description

Product Name	FILA Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC-p;IF
Species Reactivity	Human;Rat;Mouse
Specificity	FILA Polyclonal Antibody detects endogenous levels of protein.
Immunogen Description	Synthesized peptide derived from part region of human protein
Conjugates	Unconjugated
Accession No.	Swiss Prot:P20930GeneID:2312
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

IHC-p 1:50-300

Images



Immunohistochemical analysis of paraffin-embedded human liver cancer.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.