

MLKL Mouse Monoclonal Antibody(8H7)

Catalog No: #44206

Package Size: #44206-1 50ul #44206-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MLKL Mouse Monoclonal Antibody(8H7)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	8H7
Isotype	IgG1
Applications	IHC
Species Reactivity	Hu
Specificity	MLKL Mouse Monoclonal Antibody(8H7) detects endogenous levels of MLKL
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Conjugates	Unconjugated
Target Name	MLKL
Other Names	MLKL
Accession No.	Swiss-Prot:Q8NB16Gene ID:197259
Calculated MW	54kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Tonsil Tissue using MLKL Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Skin Tissue using MLKL Mouse mAb diluted at 1:200.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.